



# Hospital Mechanical Hilo Bed (Double Fowler)



## BA-3020-PN

- Double Fowler \*C/W PE Blow head & foot board
- Steel Construction
- With adjustable Height ( HI-LO Position ) Backrest & Kneerest by mean of 3 Crank Winding mechanism at foot end of bed
- Mobile on 4 x 125mm wheel Diagonally fitted w/Brakes
- Provision for IV Drip pole attachment
- Finishing of hard epoxy powder coating
- HILO Position from 420mm - 720mm ( w/o mattress )
- Mattress Size : L 1905mm x W 864 mm \* DIM : LI 2090 x W 955 x H 420-740 mm

### Key Specifications :

Hilo position from floor to mattress support top : 420-740mm

Backrest : 0° - 80°

Kneerest : 0° - 25°

Size : 211(L) x 95.5(W) cm

